

Stereo 2 – Nikon SMZ18

Quick-start Guide

Location: Room 2.136

(Updated: 31/8/2021)

Acquisition Software:

NIS Elements AR



Available Objective Lenses: 0.5x and 1.0x lens plus Magnification Range of 0.75 – 13.5x

Be mindful of what you are doing at ALL times: The microscope and the lenses are delicate and very expensive and should be treated with care.

Sample Cleaning and Preparation:


- Make sure your sample container (e.g. your Petri dish) is scrupulously clean and dry before placing it onto the microscope
- Use an appropriate sized sample holder for your sample to minimise spill risk
- Be careful not to over-fill sample containers as the microscope must be kept dry

Start-up Procedure

- Check the microscope and table are clean before you start doing anything else
- Rotate the desired objective into place and set the zoom to the lowest setting
 - o Line the objective up for the required mode of imaging – a single white dot indicates monocular while dual white dots indicate stereo mode
- Check the eyepieces to be sure that the Dioptre adjustment on both eyepieces are set to zero and that the eyepieces are a comfortable distance apart
- Pull the eyepiece/camera rod all the way out to send 100% of light to the binoculars **(1)**
- Log in to the computer using your IMB username and password
- If you are planning on doing Fluorescence, click the **Lamp ON** Icon on the desktop
 - o If the lamp won't turn on, check the control box to the left of the scope – if the yellow light is on then you will need to wait for it to turn off before starting the lamp
- Double-click the NIS-Elements icon on the desktop – this will turn on the microscope hardware, initialise the Z-Drive for the motorised focus and launch the acquisition software
 - o Make sure to let this finish completely before proceeding, failure to do so can lead to problems with Z-stack settings later in the acquisition process
- Place your sample on the stage and centre it below the objective
- Focus on your sample through the eyepieces using the focus wheel on the Nikon Remote Controller – adjust the speed with the button on the side to the rear of the focus controller

Light Path Settings

- To use the fluorescence light source, open the shutter using the puck shutter button **(2a)**
 - o Adjust the intensity of the fluorescence light using the ND selection buttons **(2b)**
 - 1 is the brightest, 32 is the dimmest
 - o To change the fluorescence filter cube, select the 1, 2 or 3 buttons on the Nikon Remote Controller **(3)** – 1 = DAPI; 2 = GFP; 3 = RFP



- To use the white light source for Brightfield imaging, turn on the lamp in the base by switching on the power at the front-left of the microscope base
 - o To adjust the intensity of the white light, rotate the control dial to the front-right of the microscope base



Saving Your Data

Save your Data to C:\UserData** and create yourself a folder if it doesn't already exist. **DO NOT** save directly to the network, only transfer your files at the end of your session.

Connecting to Network attached Storage

- o Data you acquire on this device should be backed up to a server for future retrieval and use.
- o The IMB has access to a number of servers both internally managed, as well as UQ-wide resources, think carefully where you would like to store your data.
- o Internal IMB Server address which will be UNAVAILABLE by end of 2020, include:
 - o **IMB Share**
 - Backed up
 - Group shared server
 - [\imbshare.imb.uq.edu.au](http://imbshare.imb.uq.edu.au)
 - o **IMB Transfers**
 - Not backed up
 - Deleted every 7 days
 - \transfers.imb.uq.edu.au
 - o **Objective Server**
 - Microscopy only
 - Not backed up
 - \objective.imb.uq.edu.au
- o UQ-wide resources:
 - o **UQ RDM**
 - Project based shared drive
 - Can be cached locally (\data.imb.uq.edu.au)
 - For more information see rdm.uq.edu.au

- Note this uses UQ-SSO credentials (type `uq\username`)

Basic Image Capture

- Push the eyepiece/camera rod all the way in **(1)** to send the light to the camera
- Select the required colour channel under the OC Channel tab
 - To acquire fluorescence channels, ensure the brightfield light is off and the shutter is open on the puck
 - To acquire brightfield images, close the fluorescence shutter using the puck and turn the brightfield light on
- Click Live to show a live view of the camera image
 - Adjust the resolution as needed (**Cam Rez**)
 - Change the Exposure time if too dark/bright or click Auto Exposure to have the computer adjust the exposure time automatically (**Exposure Settings**)
 - You may also adjust the display curve at the bottom of the software window to increase the image brightness or use the auto-display adjust button to do this on the fly or as a one-off snap-shot adjustment
- Click Freeze if you want to pause the live image
- Click Capture to capture the current channel, this will need to be manually saved
- Click AutoCapture and Save to capture an image of the current channel and save it automatically



Advanced Image Capture



ND Acquisition

- Set up your individual channel settings as for a basic channel setup above, these settings will be imported into the experiment acquisition setting
- Set your Experiment name at the top of the ND Acquisition window
- Ensure “Save to File” is ticked and set both your Path for the AutoSave as well as the Filename to save as
 - o While you can save as a TIFF series it is recommended to always save as a *.nd2 file format
 - o The filename will have an auto-incrementing number placed at the end of what you assign so you may capture multiple files without fear of overwriting previous ones
- To add time-lapse settings, tick the “Time” tab
 - o Add in a new time series line and set the interval and length
 - o You can add multiple time-series for complicated experimental setups
- There is no multi-position available as this system does not have a motorised stage
- To add a 3d stack tick the “Z Series” tab
 - o Set the mode of z-stack (i.e. top-bottom, centre-point or offset centre-point)
 - o Set the Top and Bottom Coordinates or the Centre point as required
 - o Set the z-step size or the slice number
- Add channels to the experimental setup by using the “Lambda” tab
 - o You must add at least one channel to your experiment to make it valid
 - o Add multiple channels if so desired
- Click **Run Now** to begin acquiring the experiment

EDF Acquisition

Extended Depth of Field acquisition allows users to acquire 3d stacks and process the data to produce a single image with all areas in focus, this is perfect for embryos and larger samples that would traditionally be hard to achieve a single image with everything in focus.

- Click the Capture EDF Image button
- Set the Z-stack settings to cover the 3d area of interest
- Untick the Use HDR button as this slows the process down massively
- Click Run to begin the automated capture
- Once finished, click the EDF Create Focussed Image to create the single plane EDF image



Shut Down Procedure

- Clean any spills within the working area and ensure the microscope is left in a clean and tidy state
- Exit the software and turn the microscope hardware off with the desktop icon
- Turn the Brightfield lamp off if you have used it
- Turn the Fluorescent lamp off via the desktop icon
- Log out of the computer – **DO NOT** shut the computer down

System Hardware Summary

Microscope Stand: Nikon SMZ18 Stereo Microscope

Objectives: 0.5x and 1.0x interchangeable lenses

Camera: Nikon Ds-Qi2 CMOS Camera

Fluorescent Filters: DAPI/GFP/RFP – Other filters are available through Imaging Facility staff.